LANEXO[®] Lab Inventory, Safety and Compliance Management System

Release Notes

LANEXO[®] System Release Notes provide a brief insight into updates on the Mobile and Web Applications.

Intended Audience

All customers and users of the LANEXO® System.

Note

If Auto-Update is disabled, please schedule an **upgrade** under **Administration menu > Version Control** on **Web app** for the latest version with new features/improvements and bug fixes.

NIGRCK

Users might have to analyze the impact for additional validation needs of the Application at their site, after the release.



The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Version History

Date	Version	Description
27-JAN-2023	HOT_2023-050	Hot fix for the Mobile Application
11-MAR-2023	R_2023-100 vADM_2.5.0 vAPP_2.1.0 vAUD_2.2.0 vINV_2.2.0 vINT_2.2.0 vLOC_2.1.0 vPRO_2.0.1 iOS Mobile App NEW	 iOS Mobile App View articles and items in an aggregated view, extended filter functionality, Consumables have been renamed to Inventory Barcode label printing from the web application for consumables and locations (to use for existing and new consumables + locations) Users are now able to configure their authentication option by logging in via access cards or email and password only Approval module is being introduced including a global LANEXO policy, users can view and work on pending approvals. Various minor improvements: Enable file deletion, allow users to set up a default value setting for life span, allow exact search in inventory, adding of additional article data fields, auto-update option for new versions Bug fixes in Administration, Audit Trail, Basket, Inventory (former Consumables), Integration, Location, System Settings, and User Profile modules.
22-APRIL-2023	Bug_2023-150	Bug fix release for the Web and Mobile Application
27-MAY-2023	Bug_2023-175	Bug fix release for the Web and Mobile Application
24-JUNE-2023	Bug_2023-180	Bug fix release for the Web Application
5-AUGUST-2023	R_2023-200 vADM_2.6.0 vAPP_2.2.0 vAUD_2.3.0 vBAS_2.1.0 vDAS_2.1.0 vEXP_2.1.0 vINV_2.3.0 vLOC_2.2.0 vPRO_2.0.2 Android Mobile App NEW	 iOS App includes now an improved Experiment module and the Inventory Count New Android App available incl. improved Experiment module New approval management is being introduced, allowing users to create custom approval policies to approval management and the inventory module Content has been updated for the English and German language (mobile and web) Module license section has been extended Bug fixes in Approval, Audit Trail, Basket, Inventory, Locations, Profile, System Settings and iOS App
25-November-2023	R_2023-250 vADM_2.6.0 vAUD_2.4.0 vBAS_2.3.0 vDAS_2.1.0 vEXP_2.2.0 vINV_2.4.0 vLOC_2.3.0 vPRO_2.0.2 Android Mobile App IOS Mobile App	 "Contact Us" button within LANEXO to enable direct Feedback to our Product Team Several Design and UI Improvements like harmonized Tables. Reason and Password for Approvers. File names for Products are now enabled to 256 Characters and Special Characters. Optimized how GHS codes are shown. Bug fixes
10-February-2024	R_2024-100 vADM_2.6.0 vAVP_2.6.0 vAUD_2.4.0 vBAS_2.3.2 vDAS_2.1.0 vEXP_2.2.2 vINV_2.5.0 vLOC_2.3.0 vPRO_2.0.2 Android Mobile App IOS Mobile App	 Inventory: Maintaining the package content during the first product registration on mobile which prevents empty items from now on Locations: Hierarchical filtering for Locations in the Inventory and Basket modules Full location path available in the application and in the file exports Substances: Filtering and showing substance data on the product aggregation table Refurbished way of maintaining substances which can now: Be added in the batch upload files for products Be linked and delinked from products Be added much easier with a new substance preview



LANEXO[®] System

Release Date: 10-February-2024

New Features

The following tables detail the new features/improvements.

Module	Feature Description			
Inventory	Be added inBe linked a	 Be added in the batch upload files for products Be linked and delinked from products 		
	Link Substance		×	
	Search by Substance Name or Ca	sNo.	Q	
		You have selected the below substance to link with the Product. You can search for new substance substance with another one. Click on "Link Substance" button to link the selected substance with		
	Substance Formula EC No Equivalent No	Hill Formula Molecular Mass Compliance		
		Ca	Incel Link Substance	

Approval

Approval policies enabled for substance linking or if a User removes a substance link.

	Approval Object	Description	Action Type	Action Behavio ur	Global Settin g	Approval Group	Approval Required	Other Conditions
Z	Substance Created during Item creation	If substance is new in the system during item creation	Notification	Non- Blocking	~	Inventory-Managers X 🗸 🖌	Users 🗸 i	Select Other Conditions
V	Substance Created manually	If substance is new in the system	Approval	Blocking	~	Inventory-Managers X 🖌	Users 🗸 i	Select Other Conditions
Y	Substance Linked	When a user links the substance to a product	Approval	Blocking	\checkmark	Inventory-Managers X 👻 👔	Users 🗸 i	Select Other Conditions
Y	Substance Unlinked	When a user removes a link from a product to a substance	Approval	Blocking	~	Inventory-Managers X 🗸 🖌	Users 🗸	Select Other Conditions
V	Substance Updated	If substance data is updated	Approval	Blocking	\checkmark	Inventory-Managers X 👻 👔	Users 👻 🚺	Select Other Conditions



Locations

Hierarchical filtering for Locations in the Inventory and Basket modules to support the User to find their Locations.

ilters				📅 Clear /	All Filters	
Product Name	SKU	Supplier		де Туре		
Type Product Name	Type SKU	Type Supplier		ose Package Type	~	
Package Content Unit Choose Package Content Unit 🗸	Package Content	Choose Location.	~			
		Clear selected ite	em			
Items per page 10 🗸 1-8 of 8 items		Building101		1 of 1 pages < 1 🗸		
Name of Product 🗠	SKU ~	Ite V Floor204	04	Status	C	
Acrylic acid (stabilised with h	8001810100	1	sition204 Container204 Ger	Reordered	0	
Art1001	1001	Loc1 1 Loc2			9	
				Rework Product		
NewProd118	1019901130	1 Merck K	GaA, Darmstadt, Ger		p	

Locations

Full location path available in the application and in the file exports to find the right Location of an Item.

> Inventory									
Product I	tems	Q Search	Exact I	Match	•	ß	驟	•	Y ¹
Filters							to c	lear All F	ilters
Product Name Type Product Name	Location Position204 < Ro	00m204 < Fl ~	Label Choose Label	~					
Items per page 10 ~ 1-1 of 1 ite	ems					1 of 1 pag	ges <	1 ~	>
Product Name 🗠	La	bel 🗠		Location	1 ∼				
Acetonitrile for DNA synthesis	S12005111AF7C			Containe	er204 < Pos	ition204	<		



Inventory

Maintaining the package content during the first product registration on mobile which prevents empty items from now on.



Note for regulated customers

With the changes on the software, the validation status of the Application needs to be reestablished by the user. Please reach out to us at <u>LANEXO@merckgroup.com</u> to receive the IQ/OQ documentation for your validation analysis.



Customer Support

In case of technical questions, please reach out to us at <u>LANEXO@merckgroup.com</u>.

Below find the customer support contact details:

Country	Email	Phone
AUSTRIA	technischerservice@merckgroup.com	+43 720380064
GERMANY	technischerservice@merckgroup.com	+49 8006271150
SWITZERLAND	technischerservice@merckgroup.com	+41 912479000
IRELAND	technicalservice@merckgroup.com	+353 1800200888
UNITED KINGDOM	technicalservice@merckgroup.com	+44 800272572
DENMARK	technicalservice@merckgroup.com	+45 82332821
FINLAND	technicalservice@merckgroup.com	+358 981710366
NORWAY	technicalservice@merckgroup.com	+47 23176070
SWEDEN	technicalservice@merckgroup.com	+46 20350511
ITALY	serviziotecnico@merckgroup.com	+39 233417340
PORTUGAL	serviciotecnico@merckgroup.com	+351 800600121
SPAIN	serviciotecnico@merckgroup.com	+34 916572043
BELGIUM	technischeservicebenelux@merckgroup.com	+32 26200260
NETHERLANDS	technischeservicebenelux@merckgroup.com	+31 207219898
FRANCE	servicetechnique@merckgroup.com	+33 800312467

SigmaAldrich.com/LANEXO

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, and LANEXO are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

